

Replacement Handset Cable

Features

- Replacement cable for the entertainment handset NIPH3-A7A
- Fitted with a teleCARE Safe Release plug
- Replacement assembly kit
- Color white



Description

Replacement cable with assembly kit for the NIPH3-A7A patient handset. This cable has a Safe Release plug and 6 connecting wires.

IMPORTANT: The Patient Handset cable should be replaced in accordance with the instructions in the Patient Handset NIPH3 Cable Replacement manual TD 92951EN. Failure to comply with those instructions could result in malfunctioning or damage to the patient handset.

Note: Ascom cannot accept responsibility for the performance of the patient handset or its ingress protection (IP) rating after the cable has been replaced.

Technical Specifications

Markets:		EU	US	UK	ANZ
Item number:	660394	•	•	•	•
Length:	2.5 meters / 98.4 inches				
Color:	White				
Number of wires:	6				
Assembly kit:	1 x adhesive cover 2 x screws 1 x lid 1 x gasket 1 x strain relief grommet 1 x cable tie				